


SYMBIOS™
GO Pump™
 RAPID RECOVERY SYSTEM



Tethered Cap

Clamp

Elastomeric Bladder

Tubing

Housing

Filter

Flow Restrictor

60 CC Syringe

Vented Cap

Male Luer Connector

Catheter

Parafix® Connector

Dressing

3 CC Syringe

Introducer Needle Sheath

Introducer & Needle

Disposible Infusion Pump Kit
 150 ml, 2 ml/hr
 Drug: _____ %
 LOT 06-100088 REF. 510042

Patient Label

Medication Label

DRUG NAME	PATIENT NAME
CONC.	DATE
VOLUME	TIME
WRAP LABEL AROUND PUMP FUSING	
INITIALS	