

Product Duration Table

| Product | Reservoir Volume | No. of Catheters | Flow Rate | Approx. Infusion Time | Standard Epidural Catheter | Fenestrated 2.5" Catheter | Fenestrated 5" Catheter | Fenestrated 7.5" Catheter | Fenestrated 10" Catheter | Units Per Case |
|-------------------------------|----------------------------|------------------|------------|-----------------------|----------------------------|---------------------------|-------------------------|---------------------------|--------------------------|----------------|
| GO Pump Single | 150ml | 1 | 2ml | 3 days | S150-E0002x1 | S150-F0252x1 | S150-F0502x1 | S150-F0752x1 | S150-F1002x1 | 5 |
| GO Pump Dual 300 | 150ml per side | 2 | 2ml | 3 days | D300-E0002x2 | D300-F0252x2 | D300-F0502x2 | D300-F0752x2 | D300-F1002x2 | 5 |
| GO Pump Dual 600 | 300ml per side | 2 | 2ml | 6 days | D600-E0002x2 | D600-F0252x2 | D600-F0502x2 | D600-F0752x2 | D600-F1002x2 | 5 |
| GO Pump 'Y'- 300 (2ml) | 300ml (150ml/ side) | 1 | 2ml | 6 days | S300-E0002x1 | S300-F0252x1 | S300-F0502x1 | S300-F0752x1 | S300-F1002x1 | 5 |
| GO Pump 'Y'- 300 (4ml) | 300ml (150ml/ side) | 1 | 4ml | 3 days | S300-E0004x1 | S300-F0254x1 | S300-F0504x1 | S300-F0754x1 | S300-F1004x1 | 5 |
| GO Pump 'Y'- 600 (6ml) | 600ml (300ml/ side) | 1 | 6ml | 4 days | S600-E0006x1 | S600-F0256x1 | S600-F0506x1 | S600-F0756x1 | S600-F1006x1 | 5 |